



MatchBox₂ series

LASERS FOR ANALYTICAL INSTRUMENTATION

- CW Lasers (Broad and Narrow Spectrum)
- Nanosecond SLM Lasers
- Accessories

ADVANTAGES

- One size
- One control interface
- One voltage (+5 VDC)
- One software
- Many wavelengths
- Low cost
- Self-contained design
- Compact and efficient
- Variety of accessories
- Dedicated for OEM

APPLICATIONS

CW BROAD SPECTRUM LASERS

- Fluorescence spectroscopy
- Scanning Microscopy
- Sorting
- Flow cytometry
- Metrology
- Optical guiding
- UV curing
- 3D printing
- Excitation


CW NARROW SPECTRUM LASERS

- Raman Spectroscopy
- Holography
- Inspection


NANOSECOND SLM LASERS

- Supercontinuum Generation
- Pulsed Laser Seeding
- Laser Induced Breakdown Spectroscopy (LIBS)
- Range Finding
- Raman Spectroscopy
- Holography


GENERAL INFORMATION



Beam profile of 1064L-11B
(far field)



Spectrum of 785L-21A SLM laser
(measurement is limited by spectrum analyzer)



Part number structure



Standard pinout

CONTROL SOFTWARE

MatchBox 2 control. 1v5. User edition.

Application settings Settings I _{LD} (max. 180mA) <input type="text" value="180"/> TEC1 temp. <input type="text" value="30"/> TEC2 temp. <input type="text" value="18"/> Optical power settings Optical power <input type="text" value="1200"/> <input checked="" type="radio"/> DAC value <input type="radio"/> mW (if calibrated)	Device functions Readings LD current 115.8mA APC TEC1 temp. 25.684 -25% TEC2 temp. N.A. 0% Body temp. 23.590 Access level 1 <input type="checkbox"/> Laser self start after power on	Device information Device found at COM5 Firmware for MatchBox II v1.6.6 Laser S/N:915302 Laser model:405L-15A 171h 8 min. 120 times
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------

Start Stop

LASER ON

MatchBox 2 series

Wave-length (nm)	Type	Output power (free space)	Output power (SM PM fiber)	Output power (MM fiber)	Wave-length tolerance +/-	Spectral linewidth FWHM (typical)	Noise (20 Hz – 20 MHz)
BROAD SPECTRUM CW LASERS							
405 nm	Diode	130 mW	40 mW	90 mW	3 nm	0.5 nm	0.5%
445 nm	Diode	50 mW	15 mW	35 mW	3 nm	0.5 nm	0.5%
488 nm	Diode	45 mW	20 mW	30 mW	3 nm	0.5 nm	0.5%
520 nm	Diode	40 mW	20 mW	30 mW	5 nm	0.5 nm	0.5%
532.1 nm	DPSS	200 mW	100 mW	160 mW	0.1 nm	0.3 nm	N/A
	DPSS	500 mW	N/A	350 mW	0.1 nm	0.3 nm	N/A
635 nm	Diode	150 mW	75 mW	120 mW	3 nm	0.5 nm	0.5%
660 nm	Diode	110 mW	45 mW	90 mW	3 nm	0.5 nm	0.5%
785 nm	Diode	150 mW	60 mW	120 mW	3 nm	0.5 nm	0.5%
830 nm	Diode	130 mW	60 mW	90 mW	10 nm	0.5 nm	0.5%
850 nm	Diode	130 mW	60 mW	90 mW	10 nm	0.5 nm	0.5%
915 nm	Diode	200 mW	80 mW	140 mW	3 nm	0.5 nm	0.5%
975 nm	Diode	200 mW	80 mW	120 mW	3 nm	0.5 nm	0.5%
980 nm	Diode	200 mW	80 mW	120 mW	3 nm	0.5 nm	0.5%
1030 nm	DPSS	500 mW	300 mW	400 mW	2 nm	2 nm	N/A
1064 nm	DPSS	500 mW	300 mW	400 mW	0.3 nm	0.5 nm	N/A
NARROW SPECTRUM (SLM) CW LASERS							
405 nm	Diode	30 mW	10 mW	20 mW	0.3 nm	<3 pm	0.5%
488 nm	Diode	30 mW	10 mW	15 mW	0.2 nm	<3 pm	0.5%
532 nm	DPSS	50 mW	25 mW	40 mW	0.3 nm	<0.2 pm	1%
632.8 nm	Diode	60 mW	30 mW	40 mW	0.1 nm	<2 pm	0.5%
635 nm	Diode	90 mW	45 mW	65 mW	0.1 nm	<2 pm	0.5%
783 nm	Diode	90 mW	45 mW	60 mW	0.1 nm	<2 pm	0.5%
785 nm	Diode	100 mW	50 mW	70 mW	0.1 nm	<2 pm	0.5%
	Diode	500 mW	N/A	350 mW	0.5 nm	<75 pm	0.5%
808 nm	Diode	100 mW	50 mW	80 mW	0.2 nm	<2 pm	0.5%
1029 nm	DPSS	400 mW	200 mW	280 mW	0.25 nm	<75 pm	1%
1064	DPSS	400 mW	200 mW	280 mW	0.3 nm	<0.2 pm	1%

Other wavelengths on request:

473 nm, 491 nm, 561 nm, 589 nm, 593 nm, 671 nm, 946 nm, 1123 nm, 1319 nm, 1342 nm.

OTHER PARAMETERS

BEAM PROPERTIES:

- Transversal mode: TEM00, except 500 mW versions of 532 nm and 785 nm
- Beam diameter at aperture (1/e²): <2 mm for diode and ~1 mm for DPSS
- Beam divergence (full angle): <2 mrad for diode and <1.5 mrad for DPSS, except 500 mW versions of 532 nm and 785 nm
- Beam pointing stability: <1 mrad/C°
- Bore sight error: +/-2 mrad (vertical), +/-3 mrad (horizontal)
- Beam quality, M2: <1.3, except multimode 500 mW versions of 532 nm and 785 nm
- Polarization ratio: better than 500:1

POWER STABILITY:

- Power stability of free-space lasers is <1 % RMS over 8 hrs
- Power stability of fiber-coupled lasers is <2 % RMS over 8 hrs
- Non-SLM DPSS lasers have significant noise peaks at above 200 kHz

MODULATION:

- Broad spectrum diode lasers can be modulated up to more than 10 MHz via TTL pin
- For SLM diode and all DPSS lasers, the TTL pin is configured for fan speed control
- DPSS lasers can be modulated at low frequencies (up to few kHz) but modulation is implemented upon request

FIBER SPECS:

- SLM fiber coupled lasers are made with FC/APC connectors
- Non-SLM lasers are made with FC/PC connectors
- Standard length of a fiber is 1 m to 1.2 m
- Polarization extinction ratio (PM fiber): better than 20 dB
- Polarization rotation (PM fiber): less than 5 degree

PHYSICAL PROPERTIES:

- Control interface type: UART serial bus, convertible to USB with standard accessories
- External power supply requirement: +5VDC, 5A for DPSS, 1.5 A for diode up to 200 mW
- Dimensions (L-W-H): 50 x 30 x 16 mm (excluding pins and output window)
- Beam height from the base: 10.4 mm (+/- 0.3 mm)
- Heatsink requirement: diode <1 °C/W, DPSS <0.5 °C/W
- Optimum heatsink temperature (non-condensing): +15...+30 °C
- Max. heatsink temperature 40 °C
- Internal temperature stabilization: TEC
- Overheat protection: Yes
- Storage temperature (non-condensing): -10 to +50 °C
- Warranty: 12 months,
- Hours limitation of 5000 hrs applies for 405, 445, 488, 515, 633, 635, 660 nm diode lasers. Operational time calculation is based on an internal EPROM counter

COMPATIBILITY:

- ROHS
- General Product Safety Directive (GPSD) 2001/95/EC
- Electromagnetic Compatibility (EMC) Directive 2004/108/EC
- IEC60825-1:2014 (compliant only using additional accessories)

Custom wavelengths and specifications are available on request

NANOSECOND PULSED LASERS



ADVANTAGES

- Same size and a physical interface as of CW MatchBox lasers.
- High pulse energy
- Single-longitudinal-mode (SLM) spectrum
- High average power
- Superb pulse-to-pulse stability
- Very low jitter

SPECIFICATION (PRELIMINARY)

Wave-length	Output Power	Pulse duration	Repetition rate	Pulse energy	Pulse-to-pulse stability (RMS)
1029 nm	500 mW	3.5 ns	20 kHz	25 μ J	<0.5 %
514.5 nm	110 mW	3.5 ns		5.5 μ J	
343 nm*	30 mW*	3.5 ns		1.5 μ J	

* Preliminary specifications

ACCESSORIES



4-port USB power supply



25 W power supply



USB power cable



Water cooled
adapter plate



Adapter to M6
taped breadboard



Air-cooled heatsink for
DPSS lasers




Protective Goggles



MatoS
Gold plated
SERS substrates



RandaS
Silver plated
SERS substrates



Contact for quote:

Phone: +370 602 24642

Email: sales@integratedoptics.com

Headquarters:

Integrated Optics, UAB

Kalvariju str. 125B, XI building

08221 Vilnius, Lithuania

Phone: +370 5 2196393

Email: info@integratedoptics.com

MatchBox™ is an unregistered trademark of Integrated Optics, UAB



www.integratedoptics.com
Made in Lithuania (EU)